

CLIC1 Antibody

Subcategory: Rabbit Polyclonal Antibody

Cat. No.: 251932

Unit: 0.1 mg

Description:

Chloride intracellular channel protein 1 (CLIC1) is a protein that can insert into membranes and form chloride ion channels. CLIC1 channel activity depends on the pH, and membrane insertion is redox-regulated and occurs only under oxidizing conditions.

Isotype: Rabbit Ig

Applications: E, WB, IHC

Species Reactivity: H


Format: Each vial contains 0.1 mg IgG in 0.1 ml (1 mg/ml) of PBS pH7.4 with 0.09% sodium azide. Antibody was purified by Protein-G affinity chromatography.

Alternate Names: Chloride intracellular channel protein 1; Nuclear chloride ion channel 27; NCC27; Chloride channel ABP; Regulatory nuclear chloride ion channel protein; hRNCC; CLIC1; NCC27

Accession No.: O00299

Antigen: KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human CLIC1.

Application Notes: E: 1:10,000-1:50,000; WB: 1:100-1:500; IHC: 1:100-1:500



The CLIC1 Antibody is used in Western blot to detect CLIC1 in human fetal kidney lysates.

Storage: Store at -20°C. Minimize freeze-thaw cycles. Product is guaranteed one year from the date of shipment.

For research use only, not for diagnostic or therapeutic procedures.