

Beta-centractin Antibody

Subcategory: Rabbit Polyclonal Antibody

Cat. No.: 251994

Unit: 0.1 mg

Description:

Beta-centractin, a member of the actin family, is a component of a multi-subunit complex involved in microtubule based vesicle motility. Beta-centractin is associated with the centrosome.

Isotype: Rabbit Ig

Applications: E, WB

Species Reactivity: H, M, R


Format: Each vial contains 0.1 mL of crude serum with 0.09% sodium azide.

Alternate Names: Beta-centractin; Actin-related protein 1B; ARP1B; ACTR1B; CTRN2

Accession No.: P42025

Antigen: KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human beta-centractin.

Application Notes: E: 1:10,000-1:50,000; WB: 1:1,000-1:2,000



The Beta-centractin Antibody is used in Western blot to detect Beta -centractin in human fetal heart lysates.

Storage: Store at -20°C. Minimize freeze-thaw cycles. Product is guaranteed one year from the date of shipment.

For research use only, not for diagnostic or therapeutic procedures.