
HSP 90-beta Antibody

Subcategory: Rabbit Polyclonal Antibody

Cat. No.: 251921

Unit: 0.1 mg

Description:

Heat shock protein HSP 90-beta (HSP 90-beta) is a molecular chaperone with ATPase activity. HSP 90-beta is a homodimer and interacts with p53/TP53. HSP 90-beta is phosphorylated upon DNA damage, probably by ATM or ATR.

Isotype: Rabbit Ig

Applications: E, WB

Species Reactivity: H

Format: Each vial contains 0.1 mg IgG in 0.1 ml (1 mg/ml) of PBS pH7.4 with 0.09% sodium azide. Antibody was purified by Protein-G affinity chromatography.

Alternate Names: Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84; HSP90AB1; HSP90B; HSPC2; HSPCB

Accession No.: P08238

Antigen: KLH-conjugated synthetic peptide encompassing a sequence within human HSP 90-beta.

Application Notes: E: 1:10,000-1:50,000; WB: 1:200-1:1,000

Storage: Store at -20°C. Minimize freeze-thaw cycles.

Product is guaranteed one year from the date of shipment.

For research use only, not for diagnostic or therapeutic procedures.