
TNFRSF19 Antibody

Subcategory: Rabbit Polyclonal Antibody

Cat. No.: 251957

Unit: 0.1 mg

Description:

Tumor necrosis factor receptor superfamily member 19 (TNFRSF19) mediates the activation of JNK and NF-kappa-B. TNFRSF19 promotes caspase-independent cell death. TNFRSF19 is highly expressed in prostate. TNFRSF19 is detected at lower levels in thymus, spleen, testis, uterus, small intestine, colon and peripheral blood leukocytes.

Isotype: Rabbit Ig

Applications: E, WB

Species Reactivity: H

Format: Each vial contains 0.1 mg IgG in 0.1 ml (1 mg/ml) of PBS pH7.4 with 0.09% sodium azide. Antibody was purified by Protein-G affinity chromatography.

Alternate Names: Tumor necrosis factor receptor superfamily member 19; Toxicity and JNK inducer; TRADE; TNFRSF19; TAJ; TROY

Accession No.: Q9NS68

Antigen: KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human TNFRSF19.

Application Notes: E: 1:10,000-1:50,000; WB: 1:100-1:500

Storage: Store at -20°C. Minimize freeze-thaw cycles.

Product is guaranteed one year from the date of shipment.

For research use only, not for diagnostic or therapeutic procedures.