

---

## **KLHL11 (1B9A) Antibody**

**Subcategory:** Mouse Monoclonal Antibody

**Cat. No.:** 253029

**Unit:** 0.1 ml

**Description:**

KLHL11 (Kelch-like protein 11) is 708 amino acid long (~80kD) secreted protein containing 5 Kelch repeats.

**Isotype:** Mouse IgG1

**Applications:** WB, IF

**Species Reactivity:** H

**Format:** Each vial contains 0.1 ml ascitic fluid with 0.03% sodium azide.

**Alternate Names:** FLJ10572; KLHL11; Kelch-like protein 11

**Accession No.:** Q9NVR0

**Antigen:** Purified recombinant fragment of human KLHL11 expressed in E. coli.

**Application Notes:** WB: 1:500-1:2000; IF: 1:200-1:1000;

**Storage:** Store at 4°C for short term use only. Store at -20°C for storage over 1 month. Product is guaranteed 6 months from the date of shipment.

For research use only, not for diagnostic or therapeutic procedures.