


# Immunoglobulins

## Antibody Solutions for Every Application

Non-immune Human IgE  
Secondary Antibodies  
Isotype Controls

### Human IgE

Cat. No. 250202 is a non-immune human IgE produced in vitro by a monoclonal hybridoma isolated from a healthy human donor. The IgE is an affinity purified fully human antibody with epsilon heavy chain and kappa light chain. The monoclonal source of IgE ensures antibody property consistency and other Ig isotype contamination free when used as standard in quantitative IgE assays. Abbiotec offers human IgE in purified, azide free and biotinylated formats.



Quality control of human IgE preparation from hybridoma (lane 2) vs. IgE sample purified from myeloma patient serum (lane 3). MW markers are shown in lane 1.

### Control and Secondary Abs

Browse our catalog for human, mouse, rabbit and rat antibodies that are developed for ELISA, western blot, immunohistochemistry and flow cytometry applications. Isotype controls, human IgE and secondary antibodies are listed below.

Cat. No.	Product Name
250202	Human IgE, 0.1 mg
250204	Human IgE, 1 mg
250203	Human IgE, azide free, 0.1 mg
250205	Human IgE, azide free, 1 mg
250206	Human IgE, biotin-conj. , 0.1 mg
250384	Anti-chicken IgY (10B10) Antibody, 0.2 mg
252240	Anti-goat IgG Antibody, HRP-conj. , 0.5 mg
250114	Anti-human IgA Antibody, 0.1 mg
250140	Anti-human IgE Antibody, 0.1 mg
250116	Anti-human IgG Antibody, 0.1 mg
250928	Anti-human IgM Antibody, 0.1 mg
250929	Anti-human IgM Antibody, PE-conj. , 2 ml
250930	Anti-human IgM Antibody, FITC-conj. , 2 ml
250987	Anti-human kappa Ig Antibody, 0.1 mg
251043	Anti-human kappa Ig Antibody, FITC-conj. , 2 ml
251099	Anti-human kappa Ig Antibody, PE-conj. , 2 ml
252241	Anti-human lambda Ig Antibody, 0.1 mg
252242	Anti-human lambda Ig Antibody, FITC-conj. , 0.1 mg
252243	Anti-human lambda Ig Antibody, PE-conj. , 2 ml
250137	Anti-mouse IgG Antibody, HRP-conj. , 1 mg
250138	Anti-mouse IgG Antibody, FITC-conj. , 1 mg
250139	Anti-mouse IgG Antibody (Fab') <sub>2</sub> , FITC-conj. , 1 mg
252237	Anti-rabbit IgG Antibody, HRP-conj. , 1 mg
252238	Anti-rabbit IgG Antibody, FITC-conj. , 1 mg
250842	Anti-rat IgA Antibody, 0.1 mg
250115	Anti-rat IgG Antibody, FITC-conj. , 0.1 mg
252239	Anti-rat IgG Antibody, HRP-conj. , 1 mg
250446	Anti-rat kappa Ig Antibody, 0.2 mg
250198	Control IgG1, mouse, 0.1 mg
250207	Control IgG1, PE-conj. , mouse, 1 ml
250199	Control IgG2a, mouse, 0.1 mg

Contact our Customer Service Dept for comprehensive listing, availability, and new product releases.



*From Biology to Discovery™*

7985 Dunbrook Rd. Ste A • San Diego, CA 92126 • USA  
1.800.854.7453 • info@abbiotec.com • www.abbiotec.com